

Product Datasheet

FAM215A antibody (orb645515)

Description	Rabbit polyclonal antibody to FAM215A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide corresponding to residues in Human FAM215A.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -25°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 12 kDa
Uniprot ID	Q9Y5M1
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:205
Expiration Date	12 months from date of receipt.