

Product Datasheet

Myc-Tag antibody (orb646361)

Description	Rabbit polyclonal antibody to Myc-Tag
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide Myc-Tag.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:2000-1:5000
Expiration Date	12 months from date of receipt.