

Product Datasheet

Actinin- alpha 1 antibody (orb646376)

Description	Rabbit polyclonal antibody to Actinin- alpha 1
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant full length Human Actinin- α 1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 103 kDa
Uniprot ID	P12814
Dilution Range	WB1:500 - 1:2000, IHC-P: -P1:50 - 1:200
Expiration Date	12 months from date of receipt.