

Product Datasheet

RAB11A antibody (orb646382)

Description	Rabbit polyclonal antibody to RAB11A
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant full length Human RAB11A.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 25 kDa
Uniprot ID	P62491
Dilution Range	WB1:500 - 1:2000, IF1:50 - 1:200
Expiration Date	12 months from date of receipt.