

Product Datasheet

EIF4G1 (phospho-S1148) antibody (orb646477)

Description	Rabbit polyclonal antibody to EIF4G1 (phospho-S1148)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide corresponding to residues in Human EIF4G1 around the phosphorylation site of S1148.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 175 kDa
Uniprot ID	Q04637
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.