

Product Datasheet

TUBB2A (Acetyl-K379) antibody (orb646485)

Description	Rabbit polyclonal antibody to TUBB2A (Acetyl-K379)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Acetylated peptide derived from human TUBB2A around the Acetylation site of K379.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 49 kDa
Uniprot ID	P68371, P04350, Q9BVA1, Q13885, Q9BUF5, Q13509, P07437
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.