

Product Datasheet

CCL25 antibody (orb646608)

Description	Rabbit polyclonal antibody to CCL25
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	KLH conjugated synthetic peptide derived from human CCL25:21-120/150
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 17 kDa
Uniprot ID	O15444
Dilution Range	WB: 1:500-1:2000, IHC-P: -P/, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.