

Product Datasheet

Histone H3 (Acetyl-K23) antibody (orb646759)

Description	Rabbit polyclonal antibody to Histone H3 (Acetyl-K23)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic peptide including acetyl-lysine 23 of Histone H3.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 17 kDa
Uniprot ID	Q71DI3, P84243, P68431
Dilution Range	WB: 1:500-1:1000, IF: 1:200
Expiration Date	12 months from date of receipt.