

Product Datasheet

Cleaved-DFF45 (D224) antibody (orb646772)

Description	Rabbit polyclonal antibody to Cleaved-DFF45 (D224)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide, corresponding Human NEIL3.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 35, 45 kDa
Uniprot ID	O00273
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.