

Product Datasheet

Collagen I antibody (orb646780)

Description	Rabbit polyclonal antibody to Collagen I
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	A synthetic peptide of human COL1A1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 139 kDa
Uniprot ID	P02452
Dilution Range	WB1:500 - 1:1000, IHC-P: -P1:50 - 1:200, IF1:50 - 1:200
Expiration Date	12 months from date of receipt.