

Product Datasheet

FBXL5 antibody (orb646944)

Description	Rabbit polyclonal antibody to FBXL5
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human FBXL5: 31-130/691
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 79 kDa
Uniprot ID	Q9UKA1
Dilution Range	WB 1:500 - 1:1000, IHC-P: -P 1:50 - 1:100
Expiration Date	12 months from date of receipt.