

Product Datasheet

GABARAPL1 antibody (orb647082)

Description	Rabbit polyclonal antibody to GABARAPL1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human GABARAPL1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 14 kDa
Uniprot ID	Q9H0R8
Dilution Range	WB1:500 - 1:2000
Expiration Date	12 months from date of receipt.