

Product Datasheet

LAMB3 antibody (orb647204)

Description	Rabbit polyclonal antibody to LAMB3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant protein of human LAMB3.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 130 kDa
Uniprot ID	Q13751
Dilution Range	IHC-P: -P1:50 - 1:100
Expiration Date	12 months from date of receipt.