

Product Datasheet

GDI1 antibody (orb647301)

Description	Rabbit polyclonal antibody to GDI1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant full length Human GDI1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 55 kDa
Uniprot ID	P31150
Dilution Range	WB1:500 - 1:2000, IHC-P: -P1:50 - 1:200
Expiration Date	12 months from date of receipt.