

Product Datasheet

TRIOBP antibody (orb647611)

Description	Rabbit polyclonal antibody to TRIOBP
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant protein of human TRIOBP
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	70kDa
Uniprot ID	Q9H2D6
Dilution Range	WB1:500 - 1:2000, IHC-P: -P1:50 - 1:200
Expiration Date	12 months from date of receipt.