

Product Datasheet

ATP6V1D antibody (orb647824)

Description	Rabbit polyclonal antibody to ATP6V1D
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant protein of human ATP6V1D
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	35kDa
Uniprot ID	Q9Y5K8
Dilution Range	WB1:500 - 1:2000, IHC-P: -P1:50 - 1:100
Expiration Date	12 months from date of receipt.