

Product Datasheet

5-Propargylamino-ddUTP (ATTO-580Q) (orb64805)

Description	5-Propargylamino-ddUTP conjugated to ATTO-580Q.
Form/Appearance	solution in water; Color: blue; pH: 7.5 ± 0.5
Concentration	0.50 mM-0.55 mM
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 587 nm, ϵ 110.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl, pH 7.5).
Formula	C ₁₂ H ₁₇ N ₃ O ₁₃ P ₃ -ATTO 580Q
Purity	≥ 95% (HPLC)
MW	Theoretical MW: 1181.20 g/mol (free acid)
SMILES	<chem>C1C[C@@H](O[C@@H]1COP(=O)(OP(=O)(OP(=O)(O)O)O)O)n1c(=O)[nH]c(=O)c(c1)C#CCN*</chem>
Expiration Date	12 months from date of receipt.