

Product Datasheet

ATP6V1A antibody (orb648057)

Description	Rabbit polyclonal antibody to ATP6V1A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human ATP6V1A
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 72 kDa
Uniprot ID	P38606
Dilution Range	WB1:500 - 1:2000
Expiration Date	12 months from date of receipt.