

Product Datasheet

PABPC1 antibody (orb648082)

Description	Rabbit polyclonal antibody to PABPC1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	A synthetic peptide of human PABPC1
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 71 kDa
Uniprot ID	P11940
Dilution Range	WB1:500 - 1:2000
Expiration Date	12 months from date of receipt.