

Product Datasheet

KLHL8 antibody (orb648086)

Description	Rabbit polyclonal antibody to KLHL8
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant protein of human KLHL8
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 69 kDa
Uniprot ID	Q9P2G9
Dilution Range	WB1:500 - 1:2000
Expiration Date	12 months from date of receipt.