

Product Datasheet

5-Propargylamino-ddUTP (5-FAM) (orb64818)

Description	5-Propargylamino-ddUTP conjugated to 5-FAM.
Form/Appearance	solution in water; Color: yellow; pH: 7.5 ± 0.5
Concentration	0.50 mM-0.55 mM
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ _{exc} 492 nm, λ _{em} 517 nm, ε 83.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
Formula	C ₃₃ H ₂₈ N ₃ O ₁₉ P ₃
Purity	≥ 95% (HPLC)
MW	Theoretical MW: 863.51 g/mol (free acid); Detected MW: 863.05 g/mol (free acid)
SMILES	<chem>C1C[C@@H](O[C@@H]1COP(=O)(OP(=O)(OP(=O)(O)O)O)O)n1c(=O)[nH]c(=O)c(c1)C#CCNC(=O)c1cc(c(cc1)c1c2ccc(=O)cc2oc2c1ccc(O)c2)C(=O)O</chem>
Expiration Date	12 months from date of receipt.