

Product Datasheet

FBLIM1 antibody (orb648345)

Description	Rabbit polyclonal antibody to FBLIM1
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human FBLIM1
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	41kDa
Uniprot ID	Q8WUP2
Dilution Range	WB1:500 - 1:2000
Expiration Date	12 months from date of receipt.