

Product Datasheet

Histone H3 (Acetyl-K27) antibody (orb648591)

Description	Rabbit polyclonal antibody to Histone H3 (Acetyl-K27)
Species/Host	Rabbit
Reactivity	Human, Mouse, Other, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC, IP, WB
Immunogen	A synthetic acetylated peptide corresponding to residues surrounding K27 of human Histone H3.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 16 kDa
Uniprot ID	P68431
Dilution Range	WB 1:1000 - 1:3000 , IHC-P: -P/, IF 1:200 - 1:500IP/, ChIP/, ChIPseq 1:20 - 1:50
Expiration Date	12 months from date of receipt.