

Product Datasheet

NEDD4L (phospho-S448) antibody (orb648601)

Description	Rabbit polyclonal antibody to NEDD4L (phospho-S448)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S448 of human NEDD4L
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	130kDa
Uniprot ID	Q96PU5
Dilution Range	WB1:500 - 1:2000
Expiration Date	12 months from date of receipt.