

Product Datasheet

Acetyl-Lysine Pan antibody (orb648608)

Description	Rabbit polyclonal antibody to Acetyl-Lysine Pan
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Acetylated BSA
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Refer to Figures
Dilution Range	WB1:500 - 1:1000IP1:20 - 1:50
Expiration Date	12 months from date of receipt.