

Product Datasheet

SLC32A1 antibody (orb648644)

Description	Rabbit polyclonal antibody to SLC32A1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant full length Human SLC32A1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~57 kDa
Uniprot ID	Q9H598
Dilution Range	WB1:1000 - 1:2000, IHC-P: -P1:50 - 1:100, IF1:50 - 1:100
Expiration Date	12 months from date of receipt.