

Product Datasheet

SKIV2L2 antibody (orb648703)

Description	Rabbit polyclonal antibody to SKIV2L2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant full length Human SKIV2L2.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~118 kDa
Uniprot ID	P42285
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200IP 1:50 - 1:200
Expiration Date	12 months from date of receipt.