

Product Datasheet

MAP1LC3A antibody (orb648705)

Description	Rabbit polyclonal antibody to MAP1LC3A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant full length Human MAP1LC3A.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~14 kDa
Uniprot ID	Q9H492 , Q9GZQ8
Dilution Range	WB: 1:500-1:2000, IHC-P: -P/, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.