

Product Datasheet

TNXB antibody (orb648728)

Description	Rabbit polyclonal antibody to TNXB
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant full length Human TNXB.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 74 kDa
Uniprot ID	P22105
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.