



Product Datasheet

NDUFV2 antibody (orb6488)

Description	NDUFV2 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human NDUFV2 (101-200/249 aa)
Target	NDUFV2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	24 kDa
Uniprot ID	P19404
Dilution Range	IHC-P=1:100-500, ICC=1:100, WB=1:500-2000, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



ICC analysis of A549 cell using NDUFV2 p...



IHC-P analysis of rat brain using NDUFV2...