

Product Datasheet

Histone H4 (sDi-Methyl R3) antibody (orb649049)

Description	Rabbit polyclonal antibody to Histone H4 (sDi-Methyl R3)
Species/Host	Rabbit
Reactivity	Human, Mouse, Other, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide conjugated to KLH derived from within residues 1-100 of Human Histone H4, sDi-methylated at R3.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 11 kDa
Uniprot ID	P62805
Dilution Range	WB1:500 - 1:2000, IHC-P: -P1:50 - 1:200, IF1:50 - 1:200
Expiration Date	12 months from date of receipt.