

Product Datasheet

5-Propargylamino-CTP (ATTO-612Q) (orb64941)

Description	5-Propargylamino-CTP conjugated to ATTO-612Q.
Form/Appearance	solution in water; pH: 7.5 ± 0.5
Concentration	1.0 mM-1.1 mM
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 615 nm, ϵ 115.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl, pH 7.5).
Formula	C ₁₂ H ₁₉ N ₄ O ₁₄ P ₃ -ATTO 612Q
Purity	≥ 95% (HPLC)
MW	Theoretical MW: 1208.22 g/mol (free acid)
SMILES	[C@@H]1([C@H]([C@@H](O[C@@H]1COP(=O)(OP(=O)(OP(=O)(O)O)O)O)n1c(=O)nc(c(c1)C#CCN*)N)O)O
Expiration Date	12 months from date of receipt.