

Product Datasheet

Ctip2 antibody (orb649692)

Description	Rabbit polyclonal antibody to Ctip2
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC, WB
Immunogen	KLH conjugated synthetic peptide derived from human Ctip2 : 331-4630/894
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 95 kDa
Uniprot ID	Q9C0K0
Dilution Range	WB:1:1, 000ICC:1:50-1:200, IHC-P: 1:50-1:200, FC:1:50-1:100
Expiration Date	12 months from date of receipt.