

Product Datasheet

Cytochrome P450 Reductase antibody [Out of stock] (orb649723)

Description	Rabbit polyclonal antibody to Cytochrome P450 Reductase
Species/Host	Rabbit
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	Recombinant protein
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
MW	77 kDa
Uniprot ID	P37040, P16435, P00388
Dilution Range	WB:1:500-1:1, 000ICC:1:50-1:200, IHC-P: 1:50-1:200, FC:1:50-1:100
Expiration Date	12 months from date of receipt.