

Product Datasheet

Estrogen Receptor beta antibody [Out of stock] (orb649830)

Description	Rabbit polyclonal antibody to Estrogen Receptor beta
Species/Host	Rabbit
Tested Applications	WB
Immunogen	recombinant protein
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
MW	59 kDa
Uniprot ID	Q92731
Dilution Range	WB:1:1, 000-1:2, 000
Expiration Date	12 months from date of receipt.