


Product Datasheet

NFkB p65 (phospho-Ser536) antibody (orb6504)

Description	NFkB p65 (phospho-Ser536) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human NFkBp65 around the phosphorylation site of Ser536 (FS(p-S)IA)
Target	phospho-NFkB p65 (Ser536)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61 kDa
Uniprot ID	Q04206
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=1µg /test
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...