

Product Datasheet

Rad51 antibody (orb650497)

Description	Rabbit polyclonal antibody to Rad51
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant protein
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	Q06609
Dilution Range	WB:1:500ICC/, IF:1:50-1:200, IHC-P: 1:50-1:200IP:1:10-1:50, FC:1:50-1:100
Expiration Date	12 months from date of receipt.