



Product Datasheet

NFKB p65 (phospho-Thr435) antibody (orb6506)


Description	NFKB p65 (phospho-Thr435) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human NFKBp65 around the phosphorylation site of Thr435 (EG(p-T)LS)
Target	phospho-NFKB p65 (Thr435)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61 kDa
Uniprot ID	Q04206
Dilution Range	WB=1:500-2000, ICC=1:100, Flow-Cyt=2ug/Test, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



Flow cytometry of Mouse spleen.Primary A...



Flow cytometry of MCF7.Primary Antibody ...



ICC analysis of Hela cell using phospho-...