

## Product Datasheet

### Tryptase alpha / beta 1 antibody [Out of stock] (orb650748)

|                                |                                                                                                          |
|--------------------------------|----------------------------------------------------------------------------------------------------------|
| <b>Description</b>             | Rabbit polyclonal antibody to<br>Tryptase alpha / beta 1                                                 |
| <b>Species/Host</b>            | Rabbit                                                                                                   |
| <b>Tested<br/>Applications</b> | IHC, IP, WB                                                                                              |
| <b>Immunogen</b>               | recombinant protein                                                                                      |
| <b>Preservatives</b>           | Rabbit IgG, 1mg/ml in PBS with<br>0.02% sodium azide, 50% glycerol,<br>pH 7.2                            |
| <b>Storage</b>                 | Store at +4°C after thawing. Aliquot<br>store at -20°C or -80°C. Avoid<br>repeated freeze / thaw cycles. |
| <b>Note</b>                    | For research use only                                                                                    |
| <b>MW</b>                      | 30 kDa                                                                                                   |
| <b>Uniprot ID</b>              | <b>P27435, Q15661, Q02844</b>                                                                            |
| <b>Dilution Range</b>          | WB:1:1, 000-1:5, 000, IHC-P: 1:50-<br>1:200                                                              |
| <b>Expiration Date</b>         | 12 months from date of receipt.                                                                          |