



Product Datasheet

NFKB p65 (phospho-Thr505) antibody (orb6508)


Description	NFKB p65 (phospho-Thr505) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human around the phosphorylation site of Thr505 (LV(p-T)G)
Target	phospho-NFKB p65 (Thr505)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61 kDa
Uniprot ID	Q04206
Dilution Range	WB=1:500-2000, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Flow cytometry of HL-60.Primary Antibody...



Western blot:Molt-4(Human) Cell;A549(Hum...



Western blot:Pancreas(Mouse);HL60(Human)...