

Product Datasheet

RPA32 (phospho T21) antibody [Out of stock] (orb650917)

Description	Rabbit polyclonal antibody to RPA32 (phospho T21)
Species/Host	Rabbit
Tested Applications	IP, WB
Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Thr21 of human RPA32/RPA2.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
MW	32 kDa
Uniprot ID	Q62193, Q63528, P15927
Dilution Range	WB:1:1, 000
Expiration Date	12 months from date of receipt.