

Product Datasheet

RSK1 p90 (Phospho-S380) antibody [Out of stock] (orb650938)

Description	Rabbit polyclonal antibody to RSK1 p90 (Phospho-S380)
Species/Host	Rabbit
Tested Applications	ICC, IP, WB
Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Ser380 of human RSK1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
MW	90 kDa
Uniprot ID	Q15418, P18653, Q63531
Dilution Range	WB:1:1, 000-1:5, 000ICC:1:50-1:200
Expiration Date	12 months from date of receipt.