

Product Datasheet

GSK3 alpha / beta (Phospho-Y216/Y279) antibody [Out of stock] (orb650980)

Description	Rabbit polyclonal antibody to GSK3 alpha / beta (Phospho-Y216/Y279)
Species/Host	Rabbit
Tested Applications	ICC, IHC, IP, WB
Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Tyr216 and 279 of human GSK3(alpha+beta).
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
MW	51 kDa
Uniprot ID	P18266, Q9WV60, P49841, Q2NL51, P49840, P18265
Dilution Range	WB:1:1, 000ICC:1:100-1:500, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.