

Product Datasheet

HSP27 (Phospho-S78) antibody [Out of stock] (orb650999)

Description	Rabbit polyclonal antibody to HSP27 (Phospho-S78)
Species/Host	Rabbit
Tested Applications	ICC, IHC, IP, WB
Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Ser78 of human Hsp27.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
MW	27 kDa
Uniprot ID	P04792
Dilution Range	WB:1:1, 000-1:2, 000ICC:1:50-1:200, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.