

Product Datasheet

TFIIB antibody (orb651211)

Description	Mouse monoclonal antibody to TFIIB
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Recombinant full length Human TFIIB.
Preservatives	Mouse IgG2b, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Predicted band size:35KDa Observed band size:37KDa
Uniprot ID	Q00403
Dilution Range	WB: 1:1000ICC: 1:50-200, IHC-P: 1:500
Expiration Date	12 months from date of receipt.