

Product Datasheet

HSP70 (N-terminus) antibody (orb651246)

Description	Mouse monoclonal antibody to HSP70 (N-terminus)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Synthetic peptide, corresponding to the N-terminus of Human HSP70.
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Predicted band size:70KDa Observed band size:70KDa
Uniprot ID	P08107
Dilution Range	WB: 1:1000ICC: 1:100-300
Expiration Date	12 months from date of receipt.