

Product Datasheet

Lck antibody (orb651286)

Description	Mouse monoclonal antibody to Lck
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant full length Human Lck.
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Predicted band size:58KDa Observed band size:58KDa
Uniprot ID	P06239
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.