

Product Datasheet

PCNA antibody (orb651294)

Description	Mouse monoclonal antibody to PCNA
Species/Host	Mouse
Reactivity	Hamster, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Recombinant full length Human PCNA.
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Predicted band size:36KDa Observed band size:36KDa
Uniprot ID	P12004
Dilution Range	WB: 1:1000ICC: 1:50-100IP: 1:50-100
Expiration Date	12 months from date of receipt.