

Product Datasheet

Strep II Tag antibody (orb651336)

Description	Mouse monoclonal antibody to Strep II Tag
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic peptide Strep II Tag.
Preservatives	Mouse IgG1, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB: 1:5000-1:10000, IF: 1:200
Expiration Date	12 months from date of receipt.