

Product Datasheet

BrdU antibody (orb651340)

Description	Mouse monoclonal antibody to BrdU
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant full length Human BrdU.
Preservatives	Mouse IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	The antibody is stable in lyophilized form if stored at -20°C or below. The reconstituted antibody can be stored for 2-3 weeks at 2-8°C. For long term storage, aliquot and store at -20°C or below. Avoid repeated freezing and thawing cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	IHC-P: 1:50-1:500, IF: 1:50-1:500, FC: 1:50-1:100
Expiration Date	12 months from date of receipt.